

For Research Use Only

SFRP1 Polyclonal antibody

Catalog Number: 28271-1-AP



Basic Information

Catalog Number: 28271-1-AP	GenBank Accession Number: BC036503	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 433 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6422	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: secreted frizzled-related protein 1	
Isotype: IgG	Calculated MW: 314 aa, 35 kDa	
Immunogen Catalog Number: AG28072	Observed MW: 30-35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human testis tissue,
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 28271-1-AP (SFRP1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.