

For Research Use Only

BPIFA2 Polyclonal antibody

Catalog Number: 28293-1-AP



Basic Information

Catalog Number:

28293-1-AP

Size:

150ul , Concentration: 800 ug/ml by Nanodrop and 367 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27401

GenBank Accession Number:

BC065726

GeneID (NCBI):

140683

UNIPROT ID:

Q96DR5

Full Name:

chromosome 20 open reading frame 70

Calculated MW:

249 aa, 27 kDa

Observed MW:

30-37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : human saliva,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human saliva tissue were subjected to SDS PAGE followed by western blot with 28293-1-AP (BPIFA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.