

For Research Use Only

RBM48 Polyclonal antibody

Catalog Number: 28307-1-AP



Basic Information

Catalog Number: 28307-1-AP	GenBank Accession Number: BC104865	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 650 ug/ml by Nanodrop;	GeneID (NCBI): 84060	Recommended Dilutions: WB 1:500-1:2000 IF/ICC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q5RL73	
Isotype: IgG	Full Name: chromosome 7 open reading frame 64	
Immunogen Catalog Number: AG27277	Observed MW: 40 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: human, mouse	WB : mouse brain tissue, mouse testis tissue IF/ICC : A431 cells,

Background Information

RBM48 (RNA Binding Motif Protein 48) is a gene encoding an RNA-binding protein, belonging to the RNA-binding motif protein family, which plays an important role in the splicing process of U12-type introns. Deletion of RBM48 leads to an increase in the retention of U12-type introns, which affects cell differentiation and proliferation. And mutations in RBM48 are associated with a variety of diseases. For example, U12 intron splicing defects are associated with myelodysplastic syndromes (MDS) and may increase the risk of acute myeloid leukemia (AML). RBM48 may serve as a prognostic marker in certain cancers, such as renal smoky cell carcinoma, renal clear cell carcinoma, and hepatocellular carcinoma.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28307-1-AP (RBM48 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using RBM48 antibody (28307-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).