

For Research Use Only

FZD1 Polyclonal antibody

Catalog Number: 28326-1-AP



Basic Information

Catalog Number:
28326-1-AP

Size:
150ul , Concentration: 133 µg/ml by
Bradford method using BSA as the
standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG28500

GenBank Accession Number:
BC051271

GeneID (NCBI):
8321

Full Name:
frizzled homolog 1 (Drosophila)

Calculated MW:
647 aa, 71 kDa

Purification Method:
Antigen affinity purification

Recommended Dilutions:
WB 1:200-1:1000

Applications

Tested Applications:
IHC, WB, ELISA

Species Specificity:
Human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 28326-1-AP (FZD1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.