

For Research Use Only

# ST6GAL2 Polyclonal antibody

Catalog Number: 28367-1-AP

1 Publications



## Basic Information

**Catalog Number:**

28367-1-AP

**Size:**

150ul, Concentration: 850 ug/ml by Nanodrop and 333 ug/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG26357

**GenBank Accession Number:**

BC008680

**GeneID (NCBI):**

84620

**UNIPROT ID:**

Q96JFO

**Full Name:**

ST6 beta-galactosamide alpha-2,6-sialyltransferase 2

**Calculated MW:**

529 aa, 60 kDa

**Observed MW:**

60 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

WB 1:300-1:800

## Applications

**Tested Applications:**

WB, ELISA

**Cited Applications:**

WB

**Species Specificity:**

Human, Mouse

**Cited Species:**

human

**Positive Controls:**

WB : HL-60 cells, RAW 264.7 cells, K-562 cells

## Notable Publications

Author	Pubmed ID	Journal	Application
Ruijia Liu	38318244	J Inflamm Res	WB

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28367-1-AP (ST6GAL2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.