

For Research Use Only

ST6GAL2 Polyclonal antibody

Catalog Number: 28367-1-AP



Basic Information

Catalog Number: 28367-1-AP	GenBank Accession Number: BC008680	Purification Method: Antigen affinity purification
Size: 150ul, Concentration: 850 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84620	Recommended Dilutions: WB 1:300-1:800
Source: Rabbit	Full Name: ST6 beta-galactosamide alpha-2,6-sialyltransferase 2	
Isotype: IgG	Calculated MW: 529 aa, 60 kDa	
Immunogen Catalog Number: AG26357	Observed MW: 60 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HL-60 cells, RAW 264.7 cells, K-562 cells
Species Specificity: Human, Mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28367-1-AP (ST6GAL2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.