

For Research Use Only

GABRB1 Polyclonal antibody

Catalog Number: 28425-1-AP



Basic Information

Catalog Number: 28425-1-AP	GenBank Accession Number: BC022449	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 µg/ml by Nanodrop and 367 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2560	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: gamma-aminobutyric acid (GABA) A receptor, beta 1	
Isotype: IgG	Calculated MW: 474 aa, 54 kDa	
Immunogen Catalog Number: AG29243	Observed MW: 50-54 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: Human, Mouse, Rat	

Background Information


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28425-1-AP (GABRB1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.