

For Research Use Only

FOXJ1 Polyclonal antibody

Catalog Number: 28443-1-AP

1 Publications



Basic Information

Catalog Number:

28443-1-AP

Size:

150ul, Concentration: 900 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28294

GenBank Accession Number:

BC046460

GeneID (NCBI):

2302

UNIPROT ID:

Q92949

Full Name:

forkhead box J1

Calculated MW:

45 kDa

Observed MW:

45-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : Caco-2 cells, HeLa cells

Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Ling Zou	36929635	Clin Respir J	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28443-1-AP (FOXJ1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.