

For Research Use Only

RSPO4 Polyclonal antibody

Catalog Number: 28559-1-AP



Basic Information

Catalog Number: 28559-1-AP	GenBank Accession Number: NM_001029871	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 750 µg/ml by Nanodrop and 380 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 343637	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: R-spondin family, member 4	
Isotype: IgG	Calculated MW: 26 kDa	
Immunogen Catalog Number: AG29308	Observed MW: ~19 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : rat brain tissue,
Species Specificity: Human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat brain tissue lysates were subjected to SDS PAGE followed by western blot with 28559-1-AP (RSPO4 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.