

For Research Use Only

WDCP Polyclonal antibody

Catalog Number: 28629-1-AP



Basic Information

Catalog Number: 28629-1-AP	GenBank Accession Number: BC035698	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 800 µg/ml by Nanodrop;	GeneID (NCBI): 80304	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: chromosome 2 open reading frame 44	
Isotype: IgG	Observed MW: 68 kDa	
Immunogen Catalog Number: AG23524		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HEK-293 cells, K-562 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28629-1-AP (C2orf44 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.