

For Research Use Only

TCTN1 Polyclonal antibody

Catalog Number: 28630-1-AP



Basic Information

Catalog Number:

28630-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30064

GenBank Accession Number:

BC062611

GeneID (NCBI):

79600

UNIPROT ID:

Q2MV58

Full Name:

tectonic family member 1

Calculated MW:

64 kDa

Observed MW:

60-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : DU 145 cells, mouse brain tissue, HCT 116 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28630-1-AP (TCTN1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.