

For Research Use Only

HRH1 Polyclonal antibody

Catalog Number: 28763-1-AP



Basic Information

Catalog Number:

28763-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30464

GenBank Accession Number:

BC060802

GeneID (NCBI):

3269

UNIPROT ID:

P35367

Full Name:

histamine receptor H1

Calculated MW:

487 aa, 56 kDa

Observed MW:

56-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-7 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells lysates were subjected to SDS PAGE followed by western blot with 28763-1-AP (HRH1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.