

For Research Use Only

# PCAF Polyclonal antibody

Catalog Number: 28770-1-AP

5 Publications



## Basic Information

Catalog Number:

28770-1-AP

Size:

150ul, Concentration: 700 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30238

GenBank Accession Number:

BC060823

GeneID (NCBI):

8850

UNIPROT ID:

Q92831

Full Name:

K(lysine) acetyltransferase 2B

Calculated MW:

93 kDa

Observed MW:

95-110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

Human

Cited Species:

human, mouse

Positive Controls:

WB : A431 cells, HeLa cells

## Notable Publications

Author	Pubmed ID	Journal	Application
Lipeng Pei	39901144	Mol Med	WB
Enyuan Zhang	39012622	J Proteome Res	WB
Yihong Li	38821061	Cancer Cell	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1%BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28770-1-AP (KAT2B antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.