

For Research Use Only

DHX30 Polyclonal antibody

Catalog Number: 28814-1-AP



Basic Information

Catalog Number:

28814-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop and 300 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29885

GenBank Accession Number:

BC014237

GeneID (NCBI):

22907

UNIPROT ID:

Q7L2E3

Full Name:

DEAH (Asp-Glu-Ala-His) box polypeptide 30

Calculated MW:

501 aa, 56 kDa

Observed MW:

134 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293T cells, HuH-7 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28814-1-AP (DHX30 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.