

For Research Use Only

DHX30 Polyclonal antibody

Catalog Number: 28814-1-AP



Basic Information

Catalog Number: 28814-1-AP	GenBank Accession Number: BC014237	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 22907	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: DEAH (Asp-Glu-Ala-His) box polypeptide 30	
Isotype: IgG	Calculated MW: 501 aa, 56 kDa	
Immunogen Catalog Number: AG29885	Observed MW: 134 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293T cells, HuH-7 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28814-1-AP (DHX30 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.