

For Research Use Only

mouse anti Halo Monoclonal antibody (28A8)



www.ptglab.com

Catalog Number: 28a8

4 Publications

Basic Information

Catalog Number:
28a8

Clone No.:
28A8

Applications:
WB

Host/IsoType:
Mouse/IgG1

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2827565

Purification Method:
Protein A purified

Description

Mouse monoclonal antibody [28A8] to Halo-tag

Specificity/Target

Halo-tag, derived from the bacterial haloalkane dehalogenase enzyme from *Rhodococcus rhodochrous*

Recommended Dilution

Western Blot: 1:500

Storage

Storage:

Shipped on blue ice. Store at -20°C. Stable for one year.

Storage Buffer:


PBS, 50% glycerol, 0.09% sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Western blot analysis of cell extract from HEK293T cells transiently expressing Actin-Chromobody-EGFP-Halo-tag and from untransfected cells.



Sensitivity tested with purified recombinant Halo-tag protein. Primary antibody: anti-Halo 1:500. Secondary antibody: anti-mouse Alexa Fluor[®] 647 1:1,000.