

For Research Use Only

DKK2 Polyclonal antibody

Catalog Number: 29023-1-AP

1 Publications



Basic Information

Catalog Number:

29023-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30596

GenBank Accession Number:

BC075078

GeneID (NCBI):

27123

UNIPROT ID:

Q9UBU2

Full Name:

dickkopf homolog 2 (Xenopus laevis)

Calculated MW:

259 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, mouse, rat

Cited Species:

rat

Positive Controls:

WB : mouse heart tissue, rat heart tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Yixuan Wang	39667489	Life Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29023-1-AP (DKK2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.