

For Research Use Only

DHX58 Polyclonal antibody

Catalog Number: 29114-1-AP



Basic Information

Catalog Number: 29114-1-AP	GenBank Accession Number: BC014949	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 240 µg/ml by Nanodrop;	GeneID (NCBI): 79132	
Source: Rabbit	Full Name: DEXH (Asp-Glu-X-His) box polypeptide 58	
Isotype: IgG	Calculated MW: 678 aa, 77 kDa	
Immunogen Catalog Number: AG30760		

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data