

For Research Use Only

SEC22A Polyclonal antibody, PBS Only

Catalog Number: 29179-1-PBS



Basic Information

Catalog Number:

29179-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30293

GenBank Accession Number:

BC007122

GeneID (NCBI):

26984

UNIPROT ID:

Q961W7

Full Name:

SEC22 vesicle trafficking protein
homolog A (*S. cerevisiae*)

Calculated MW:

35 kDa

Observed MW:

30-35 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29179-1-AP (SEC22A antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 29179-1-PBS in a different storage buffer formulation.