

For Research Use Only

ND3 Polyclonal antibody

Catalog Number: 29341-1-AP



Basic Information

Catalog Number: 29341-1-AP	GenBank Accession Number: YP_003024033	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 170 µg/ml by Nanodrop;	GeneID (NCBI): 4537	Recommended Dilutions: IF 1:200-1:800
Source: Rabbit	Full Name: NADH dehydrogenase, subunit 3 (complex I)	
Isotype: IgG	Calculated MW: 13 kDa	
Immunogen Catalog Number: AG29115		

Applications

Tested Applications: IF, ELISA	Positive Controls: IF : HEK-293 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using ND3 antibody (29341-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).