

For Research Use Only

WDR47 Polyclonal antibody

Catalog Number: 29364-1-AP



Basic Information

Catalog Number:

29364-1-AP

Size:

150ul , Concentration: 900 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30853

GenBank Accession Number:

BC039254

GeneID (NCBI):

22911

UNIPROT ID:

O94967

Full Name:

WD repeat domain 47

Calculated MW:

919 aa, 102 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, K-562 cells

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29364-1-AP (WDR47 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.