

For Research Use Only

WDR47 Polyclonal antibody

Catalog Number: 29364-1-AP



Basic Information

Catalog Number: 29364-1-AP	GenBank Accession Number: BC039254	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop;	GeneID (NCBI): 22911	Recommended Dilutions: WB 1:2000-1:16000
Source: Rabbit	Full Name: WD repeat domain 47	
Isotype: IgG	Calculated MW: 919 aa, 102 kDa	
Immunogen Catalog Number: AG30853	Observed MW: 100 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, K-562 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29364-1-AP (WDR47 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.