

For Research Use Only

SOX11 Polyclonal antibody

Catalog Number: 29395-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

29395-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30020

GenBank Accession Number:

NM_003108

GeneID (NCBI):

6664

UNIPROT ID:

P35716

Full Name:

SRY (sex determining region Y)-box 11

Calculated MW:

47 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : fetal human brain tissue, Y79 cells

Notable Publications

Author	Pubmed ID	Journal	Application
Shaoqin Tu	36910334	Stem Cells Int	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29395-1-AP (SOX11 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.