

For Research Use Only

SOX11 Polyclonal antibody

Catalog Number: 29395-1-AP **1 Publications**



Basic Information

Catalog Number: 29395-1-AP	GenBank Accession Number: NM_003108	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6664	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: SRY (sex determining region Y)-box 11	
Isotype: IgG	Calculated MW: 47 kDa	
Immunogen Catalog Number: AG30020	Observed MW: 60 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : fetal human brain tissue, Y79 cells
Cited Applications: WB	
Species Specificity: Human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shaoqin Tu	36910334	Stem Cells Int	WB

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29395-1-AP (SOX11 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.