

For Research Use Only

SCNN1A Polyclonal antibody

Catalog Number: 29506-1-AP



Basic Information

Catalog Number:

29506-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31586

GenBank Accession Number:

NM_001038

GeneID (NCBI):

6337

UNIPROT ID:

P37088

Full Name:

sodium channel, nonvoltage-gated 1

alpha

Calculated MW:

76KD

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

