

For Research Use Only

RASGRP4 Polyclonal antibody, PBS Only

Catalog Number: 29510-1-PBS



Basic Information

Catalog Number: 29510-1-PBS	GenBank Accession Number: BC146669	Purification Method: Antigen affinity purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 115727	
Source: Rabbit	UNIPROT ID: Q8TDF6	
Isotype: IgG	Full Name: RAS guanyl releasing protein 4	
Immunogen Catalog Number: AG31219	Calculated MW: 673 aa, 75 kDa	
	Observed MW: 54-75 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

RASGRP4 is a member of the RAS guanine nucleotide-releasing protein (RASGRP) family, which are crucial guanine nucleotide exchange factors (GEFs). The primary function of this family is to act as molecular switches that activate RAS signaling pathways by promoting the exchange of GDP for GTP on small GTPases. Specifically, RASGRP4 is highly and preferentially expressed in mast cells. It plays a vital role in the development, maturation, and functional activation of these immune cells by regulating key signaling cascades, such as the MAPK pathway.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29510-1-AP (RASGRP4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 29510-1-PBS in a different storage buffer formulation.