

For Research Use Only

# SH3BP4 Polyclonal antibody

Catalog Number: 29623-1-AP



## Basic Information

Catalog Number:

29623-1-AP

Size:

150ul, Concentration: 600 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31221

GenBank Accession Number:

BC057396

GeneID (NCBI):

23677

UNIPROT ID:

Q9POV3

Full Name:

SH3-domain binding protein 4

Calculated MW:

963 aa, 108 kDa

Observed MW:

110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : HepG2 cells, MCF-7 cells, mouse lung tissue, rat lung tissue

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29623-1-AP (SH3BP4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.