

For Research Use Only

UBAP2 Polyclonal antibody

Catalog Number: 29628-1-AP



Basic Information

Catalog Number: 29628-1-AP	GenBank Accession Number: NM_018449	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 ug/ml by Nanodrop;	GeneID (NCBI): 55833	Recommended Dilutions: WB 1:5000-1:50000 IF/ICC 1:200-1:800
Source: Rabbit	UNIPROT ID: Q5T6F2	
Isotype: IgG	Full Name: ubiquitin associated protein 2	
Immunogen Catalog Number: AG30572	Calculated MW: 117 kDa	
	Observed MW: 117 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: human	WB : HeLa cells, IF/ICC : U2OS cells,

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29628-1-AP (UBAP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed U2OS cells using UBAP2 antibody (29628-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).