

For Research Use Only

UBAP2 Polyclonal antibody

Catalog Number: 29628-1-AP



Basic Information

Catalog Number: 29628-1-AP	GenBank Accession Number: NM_018449	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop;	GeneID (NCBI): 55833	Recommended Dilutions: WB 1:5000-1:50000
Source: Rabbit	Full Name: ubiquitin associated protein 2	
Isotype: IgG	Calculated MW: 117 kDa	
Immunogen Catalog Number: AG30572	Observed MW: 117 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29628-1-AP (UBAP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.