

For Research Use Only

C12orf45 Polyclonal antibody, PBS Only

Catalog Number: 29660-1-PBS



Basic Information

Catalog Number: 29660-1-PBS	GenBank Accession Number: BC032326	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 121053	
Source: Rabbit	UNIPROT ID: Q8N5J9	
Isotype: IgG	Full Name: chromosome 12 open reading frame 45	
Immunogen Catalog Number: AG31258	Calculated MW: 188 aa, 20 kDa	
	Observed MW: 25-30 kDa	

Applications

Tested Applications:
IP, Indirect ELISA

Species Specificity:
human

Background Information

C12ORF45, also known as NOPCHAP1, is an assembly factor for the NOP-protein chaperone 1, which is an intrinsically disordered protein in free state. NOPCHAP1 plays a bridging role in the assembly of C/D snoRNPs, facilitating the selection of NOP58 and the assembly of box C/D snoRNPs.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-C12orf45 (IP:29660-1-AP, 4ug; Detection:29660-1-AP 1:400) with HEK-293T cells lysate 1320 ug. This data was developed using the same antibody clone with 29660-1-PBS in a different storage buffer formulation.