

For Research Use Only

CDT1 Polyclonal antibody

Catalog Number: 29766-1-AP



Basic Information

Catalog Number: 29766-1-AP	GenBank Accession Number: NM_030928	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 µg/ml by Nanodrop;	GeneID (NCBI): 81620	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: chromatin licensing and DNA replication factor 1	
Isotype: IgG	Calculated MW: 60KD	
Immunogen Catalog Number: AG31758	Observed MW: 65-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MCF-7 cells, HeLa cells, HepG2 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29766-1-AP (CDT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.