

For Research Use Only

CDT1 Polyclonal antibody

Catalog Number: 29766-1-AP



Basic Information

Catalog Number:

29766-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31758

GenBank Accession Number:

NM_030928

GeneID (NCBI):

81620

UNIPROT ID:

Q9H211

Full Name:

chromatin licensing and DNA replication factor 1

Calculated MW:

60KD

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : MCF-7 cells, HeLa cells, HepG2 cells

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29766-1-AP (CDT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.