

For Research Use Only

SPINK5 Polyclonal antibody

Catalog Number: 29808-1-AP

1 Publications



Basic Information

Catalog Number: 29808-1-AP	GenBank Accession Number: NM_001127698	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 11005	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9NQ38	
Isotype: IgG	Full Name: serine peptidase inhibitor, Kazal type 5	
Immunogen Catalog Number: AG30432	Calculated MW: 121 kDa	
	Observed MW: 120-130 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COLO 320 cells,
Species Specificity: Human	

Notable Publications

Author	Pubmed ID	Journal	Application
Yaning Guan	41080545	Front Immunol	

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 29808-1-AP (SPINK5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.