

For Research Use Only

WDR4 Polyclonal antibody

Catalog Number:30062-1-AP



Basic Information

Catalog Number: 30062-1-AP	GenBank Accession Number: BC001074	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 460 ug/ml by Nanodrop;	GeneID (NCBI): 10785	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	UNIPROT ID: P57081	
Isotype: IgG	Full Name: WD repeat domain 4	
Immunogen Catalog Number: AG30580	Calculated MW: 30 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse kidney tissue, rat kidney tissue
Species Specificity: human, mouse, rat	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30062-1-AP (WDR4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.