

For Research Use Only

IL-1 alpha Polyclonal antibody

Catalog Number:30216-1-AP



Basic Information

Catalog Number: 30216-1-AP	GenBank Accession Number: NM_010554	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 ug/ml by Nanodrop;	GeneID (NCBI): 16175	Recommended Dilutions: WB: 1:500-1:1000
Source: Rabbit	UNIPROT ID: P01582	
Isotype: IgG	Full Name: interleukin 1 alpha	
Immunogen Catalog Number: EG0263	Calculated MW: 31kd	
	Observed MW: 31 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : LPS treated RAW 264.7 cells,
Species Specificity: mouse	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



LPS treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 30216-1-AP (IL1-alpha antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.