

For Research Use Only

SLC52A2 Polyclonal antibody, PBS Only

Catalog Number: 30396-1-PBS



Basic Information

Catalog Number: 30396-1-PBS	GenBank Accession Number: BC002917	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 79581	
Source: Rabbit	UNIPROT ID: Q9HAB3	
Isotype: IgG	Full Name: Solute carrier family 52, riboflavin transporter, member 2	
Immunogen Catalog Number: AG33252	Calculated MW: 46 kDa	

Applications

Tested Applications:
IHC, Indirect ELISA

Species Specificity:
human

Background Information

SLC52A2 (Also known as RFVT2, RFT3) is a membrane protein that belongs to the riboflavin transporter family. SLC52A2 is a transmembrane protein that mediates cellular uptake of riboflavin. Humans are unable to synthesize vitamin B2/riboflavin and must obtain it via intestinal absorption (PubMed:20463145). Mutations in SLC52A2 have been associated with Brown-Vialetto-Van Laere syndrome 2 (BVVLS2) - an autosomal recessive progressive neurologic disorder characterized by deafness, bulbar dysfunction, and axial and limb hypotonia (PMID: 22864630, 23243084).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 30396-1-AP (SLC52A2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 30396-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 30396-1-AP (SLC52A2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 30396-1-PBS in a different storage buffer formulation.