

For Research Use Only

INCA1 Polyclonal antibody

Catalog Number: 30441-1-AP



Basic Information

Catalog Number:

30441-1-AP

Size:

150ul , Concentration: 550 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27997

GenBank Accession Number:

BC119781

GeneID (NCBI):

388324

UNIPROT ID:

Q0VD86

Full Name:

inhibitor of CDK interacting with cyclin A1

Calculated MW:

236 aa, 27 kDa

Observed MW:

26-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30441-1-AP (INCA1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.