

For Research Use Only

SEC23A Polyclonal antibody

Catalog Number:30465-1-AP



Basic Information

Catalog Number: 30465-1-AP	GenBank Accession Number: BC036649	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 ug/ml by Nanodrop;	GeneID (NCBI): 10484	Recommended Dilutions: WB: 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q15436	
Isotype: IgG	Full Name: Sec23 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG33185	Calculated MW: 86 kDa	
	Observed MW: 86 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COS-7 cells, NIH/3T3 cells, PC-12 cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30465-1-AP (SEC23A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.