

For Research Use Only

MTMR9 Polyclonal antibody

Catalog Number: 30468-1-AP



Basic Information

Catalog Number:

30468-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33501

GenBank Accession Number:

BC022003

GeneID (NCBI):

66036

UNIPROT ID:

Q96QG7

Full Name:

myotubularin related protein 9

Calculated MW:

63 kDa

Observed MW:

60-63 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HepG2 cells, mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30468-1-AP (MTMR9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.