

For Research Use Only

DUSP10 Polyclonal antibody

Catalog Number: 30493-1-AP



Basic Information

Catalog Number: 30493-1-AP	GenBank Accession Number: BC031405	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 ug/ml by Nanodrop;	GeneID (NCBI): 11221	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9Y6W6	
Isotype: IgG	Full Name: dual specificity phosphatase 10	
Immunogen Catalog Number: AG32718	Observed MW: 53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, mouse heart tissue, rat brain tissue
Species Specificity: Human, Mouse, Rat	

Background Information

DUSP10, also known as MKP5, is a member of the MKPs subfamily involved in cell proliferation, differentiation, and migration. DUSP10 is a dual specificity phosphatase capable of dephosphorylating both pTyr and pThr residues of activated MAPKs with different affinities. It has been shown that isoforms of the p38 and JNK subfamilies are selectively and more effectively dephosphorylated than ERK subfamily members. DUSP10 has been proposed as a target for therapeutic treatment of atherosclerosis, since DUSP10 inhibition reduces NF- κ B-induced TNF- α expression and increases TGF- β 1 levels in this disease (PMID: 30939861).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30493-1-AP (DUSP10 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.