

For Research Use Only

SRRM2 Polyclonal antibody

Catalog Number: 30741-1-AP



Basic Information

Catalog Number: 30741-1-AP	GenBank Accession Number: NM_016333.3	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop;	GeneID (NCBI): 23524	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: serine/arginine repetitive matrix 2	
Isotype: IgG	Calculated MW: 300 kDa	
Immunogen Catalog Number: AG27270		

Applications

Tested Applications: IF, ELISA	Positive Controls: IF : U2OS cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using SRRM2 antibody (30741-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).