

For Research Use Only

PSD Polyclonal antibody

Catalog Number: 30860-1-AP



Basic Information

Catalog Number:

30860-1-AP

Size:

150ul , Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34130

GenBank Accession Number:

BC142643

GeneID (NCBI):

5662

UNIPROT ID:

A5PKW4

Full Name:

pleckstrin and Sec7 domain containing

Observed MW:

130-140kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse , rat, human

Positive Controls:

WB : rat brain tissue, mouse brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30860-1-AP (PSD antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.