

For Research Use Only

MPC1 Polyclonal antibody, PBS Only

Catalog Number:30868-1-PBS



Basic Information

Catalog Number: 30868-1-PBS	GenBank Accession Number: BC000810	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51660	
Source: Rabbit	UNIPROT ID: Q9Y5U8	
Isotype: IgG	Full Name: brain protein 44-like	
Immunogen Catalog Number: AG28826	Calculated MW: 12 kDa	
	Observed MW: 12 kDa, 20-24 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

MPC1, also known as BRP44L. BRP44L is predicted to be located in Mitochondrion inner membrane. The protein forms a heterodimer with MPC2, and the heterodimer is a more stable and dominant form (PMID:26253029, PMID:32403431). Mpc1 and Mpc2 associate to form an ~150-kilodalton complex in the inner mitochondrial membrane. Yeast and Drosophila mutants lacking MPC1 display impaired pyruvate metabolism, with an accumulation of upstream metabolites and a depletion of tricarboxylic acid cycle intermediates. Loss of yeast Mpc1 results in defective mitochondrial pyruvate uptake, and silencing of MPC1 or MPC2 in mammalian cells impairs pyruvate oxidation. A point mutation in MPC1 provides resistance to a known inhibitor of the mitochondrial pyruvate carrier (PMID: 22628558).The molecular weight of MPC1 is 12 kDa.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30868-1-AP (BRP44L antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30868-1-PBS in a different storage buffer formulation.