

For Research Use Only

SHH Polyclonal antibody

Catalog Number: 30951-1-AP



Basic Information

Catalog Number:

30951-1-AP

Size:

150ul , Concentration: 300 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009170.3

GeneID (NCBI):

20423

UNIPROT ID:

Q62226

Full Name:

sonic hedgehog

Calculated MW:

48KD

Observed MW:

19-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : mouse brain tissue, mouse lung tissue, mouse stomach tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30951-1-AP (SHH antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.