

For Research Use Only

AMAC1L2 Polyclonal antibody, PBS Only

Catalog Number:30979-1-PBS



Basic Information

Catalog Number: 30979-1-PBS	GenBank Accession Number: BC131696	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 83650	
Source: Rabbit	UNIPROT ID: Q96KT7	
Isotype: IgG	Full Name: acyl-malonyl condensing enzyme 1-like 2	
Immunogen Catalog Number: AG34353	Calculated MW: 388 aa, 35 kDa	
	Observed MW: 34 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

AMAC1L2, also known as SLC35G5, belongs to the transporter protein family 35 member G5. This protein has an N-terminal α -helix, a central β -side and a C-terminal structural domain at the cell membrane. The role of SLC35G5 in the tumor microenvironment as a drug target and also as a biomarker for diagnosis and treatment monitoring.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30979-1-AP (AMAC 1L2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30979-1-PBS in a different storage buffer formulation.