

For Research Use Only

VEGFB Polyclonal antibody

Catalog Number: 30980-1-AP



Basic Information

Catalog Number:

30980-1-AP

Size:

150ul , Concentration: 900 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34269

GenBank Accession Number:

BC008818

GeneID (NCBI):

7423

UNIPROT ID:

P49765

Full Name:

vascular endothelial growth factor B

Calculated MW:

22 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein Eg0123 were subjected to SDS PAGE followed by western blot with 30980-1-AP (VEGFB antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.