

For Research Use Only

PTF1A Polyclonal antibody

Catalog Number: 30987-1-AP



Basic Information

Catalog Number:

30987-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34321

GenBank Accession Number:

NM_178161

GeneID (NCBI):

256297

UNIPROT ID:

Q7RTS3

Full Name:

pancreas specific transcription factor, 1a

Calculated MW:

35 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat, human

Positive Controls:

WB : mouse pancreas tissue, rat pancreas

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30987-1-AP (PTF1A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.