

For Research Use Only

MFRP Polyclonal antibody

Catalog Number:31017-1-AP



Basic Information

Catalog Number:

31017-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34302

GenBank Accession Number:

BC160191

GeneID (NCBI):

83552

UNIPROT ID:

Q9BY79

Full Name:

membrane frizzled-related protein

Calculated MW:

62 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : fetal human brain tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31017-1-AP (MFRP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.