

For Research Use Only

# Collagen Type XV Polyclonal antibody, PBS Only

Catalog Number:31056-1-PBS



## Basic Information

Catalog Number:

31056-1-PBS

Size:

100ug , Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34644

GenBank Accession Number:

BC157094

GeneID (NCBI):

1306

UNIPROT ID:

P39059

Full Name:

collagen, type XV, alpha 1

Observed MW:

250 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31056-1-AP (Collagen Type XV antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 31056-1-PBS in a different storage buffer formulation.