

For Research Use Only

C18orf34 Polyclonal antibody

Catalog Number: 31148-1-AP



Basic Information

Catalog Number: 31148-1-AP	GenBank Accession Number: BC130535	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 500 ug/ml by Nanodrop;	GeneID (NCBI): 374864	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q5BJE1	
Isotype: IgG	Full Name: chromosome 18 open reading frame 34	
Immunogen Catalog Number: AG33847	Observed MW: 95-102 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse spleen tissue, rat ovary tissue, rat spleen tissue
Species Specificity: Human, mouse, rat	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31148-1-AP (C18orf34 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.